ENTRALTM

SONDE NASO-GASTRIQUE D'ALIMENTATION ENTERALE
NASOGASTRIC ENTERAL FEEDING TUBE
SONDA D'ALIMENTAZIONE ENTERALE NASO GASTRICA.
SONDA NASOGÁSTRICA PARA ALIMENTACIÓN ENTERAL
ENTERALE NASALE MAGENSONDE

Description :

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Description:

The ENTRAL feeding tubes are designed for nasogastric or nasointestinal enteral feeding. The tubes are made of biocompatible polyurethane and are radiopac or have a line which can be localised by X-rays. Certain references have a stainless steel stylet to make the insertion easier. The weighted tubes have a distal end of tungsten.
TUBE EQUIPPED WITH A SAFETY CONNECTOR LUER LOCK MALE in line with EN 1615 European Standard.

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INSTRUCTIONS FOR USE:

Warning: The following procedure is intended for reference only, particular situations may require modifications.

1 Explain procedure to patient.

2 Place the patient in sitting position, leaning slightly backwards.

3 Measure tube for placement. Extend tube from the patient's earlobe to the bridge of the nose and from there to the tip of the xiphoid process. Reference marks are located every 10 cm from the distal end of the tube.

4 ENTRAL tubes are prelubricated and do not need any particular preparation before insertion. It is possible to improve lubrication by applying some silicone gel on the tube.

5 Insert the tube following hospital procedure. Swallowing makes insertion easier (with or without liquid absorption).

6 Before removing the stylet the position of the tube must be checked

by:

- X-ray control- fiber optics- aspiration of stomach juice - stethoscope during air injection

7 Remove stylet.

8 Connect tube to administration set and start feeding.

CAUTION

Single use product, do not reuse. Do not insert the stylet into the tube without having removed the tube entirely from the patient. Use only ENTRAL stylets. The tube must be rinsed as soon as feeding is interrupted. Use only syringes with volumes of at least 50 cc's to check the correct position and for the rinsing of the tube. Replace tube after 3 to 4 weeks.

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